

# SAFE-MED

Pressure Distributing Products  
www.safemed.dk

To order cover in Soft-Cell™ white colour or in Comfor™ PUR light grey colour.  
Cotton velour white colour is also available for Ear Pillow.

No. 115/HMI-no. 54575

Ear Pillow for pressure relief of the ear lying on the side or as pillow for pressure relief of the back of the head lying on the back.  
40x32x6 cm, hole Ø 10 cm



No. 106.1/HMI-no. 54420, 45x30x3/11 cm

No. 106.2/HMI-no. 89682, 45x40x3/11 cm


No. 106.3/HMI-no. 120870, 80x30x3/12 cm

No. 106.3518/HMI-no. 129051, 35x18x3/8 cm


Support wedge for 30 degree side positioning.  
Seat wedge for under thighs to avoid shearing forces while sitting in bed.



SAFE-MED  
Trykfordelende Produkter  
www.safemed.dk



## Pressure relief cushions



No. 104YP.2515/HMI-no. 73184

1 set Soft-Lift cushions for

"Micro-stimulation"

2 pcs. cushions 25x15x6 cm



Positioning with Soft-Lift cushions

"Micro-stimulation"

To use at points under the body for pressure distribution/pressure relief and thus pain relief.

Can be used without cover or with

Soft-Cell white cover no. 104TEN.2515