


COMFOR™ OP Cover

- black color

- Cover for cushions & mattress to the OP-theater, 210g/m²
- Polyester 125g, PU coating 85 g
Free from PVC, latex and phthalates
Biocompatible to ISO 10993-5
- Tested as medical device
- Strong, elastic, breathable
- Antimicrobial
- X-ray penetrable and it is harmless by occasional contact
- Flame test BS 7175, Crib 5
Ignition and match flame test EN 597-1 & 2
- Machine wash up to 95°C / test 50 times
- Tumble dry up to 60°C / test 50 times
- Wipe off with soapy water, alcohol solution or chlorine disinfection or similar with active chlorine 1000-1200 ppm
- 200 cm water-column(waterproof)
- CE-marked Medical Device MDR regulation (EU) 2017/745
- Permanently anti-static BS 2050:≤ 10⁶ Ohm
- The anti-static property is best in black color, which is why all covers are black. Color mark on the cover is possible and can be ordered as desired with a Ø 2 cm mark in green, orange, yellow, blue or red.



**Comfor OP covers have concealed plastic zipper and glued seams.
Can be supplied with anti-slip backside**

