

**Choose the type of cover for SAFE Med overlay mattresses and fullbody mattresses in the desired size.**

Soft-Cell or Comfor PUR Incontinence cover as full cover with zipper or as a cover with open back side.

See product details and prices and state the desired cover type in connection with ordering.

- Soft-Cell™ waterproofcover
- With fabric surface like SILK
- Fabric surface with 30g PU waterproofbackside
- Strong, elastic, anti-static, breathableand absorbs more moisture than cotton
- Ignition test EN 597-1:1995 & ISO 12952-1:2011
- Washable in machine 60°C-95°C
- By 60°C – at least 100 times
- By 95°C – at least 20 times
- Tumble dry up to 70C
- Breathability1000 g/m2 ·24h
- Öko TEX Class IProduct
- Comfor™ PUR Incontinence cover
- PU waterproofsurface on polyester
- Strong, elastic, breathable
- Antimicrobiel
- Flametest BS 7175, Crib 5
- Ignition & match flame test EN 597-1 & 2
- Washing in machine up to 95°C
- Tumble dry up to 60C
- Test 50 times wash and tumble drying
- Breathability500 g/m2 ·24h
- Wipe off with soapy water, alcohol solution orWet wipe similar to mild chlorine 1000-1200 ppm
- Biocompatible to ISO 109935

**\*All covers are FREE from PVC, latex and phthalates**

**Prices on Soft-Cell covers are ex. VAT:**

**Soft-Cell full cover with concealed zipper on 3 sides and anti-slip back side with height 3,4 or 6 cm.**

80,85,90,100 and 120 cm width DKK. 868,-  
140 cm width DKK. 923,-

**Soft-Cell cover with open back side with height 3,4 or 6 cm**

80,85 width DKK. 381,-  
90 cm width DKK. 408,-  
100 cm width DKK. 422,-  
120, 140 cm width DKK. 526,-  
180 cm width DKK. 876,-

**Soft-Cell cover with open back side with height 20 or 30 cm**

90 cm width DKK. 754,-  
140 cm width DKK. 1254,-

**Soft-Cell sheet**

220 x 140 cm DKK. 526,-

Standard colour is white.

Colour grey or black can be chosen in Soft-Cell with a surcharge of DKK. 50,-

**Prices on Comfor PUR incontinence covers are ex. VAT:**

**Comfor PUR incontinence cover with concealed zipper on 3 sides and anti-slip backside with height 3,4 or 6 cm.**

80,85,90,100 and 120 cm width DKK. 868,-  
140 cm width DKK. 923,-

**Comfor PUR incontinence cover with open back side with height 3,4 or 6 cm**

80,85 width DKK. 622,-  
90 cm width DKK. 622,-  
100 cm width DKK. 622,-  
120, 140 cm width DKK. 754,-  
180 cm width DKK. 1104,-

**Comfor PUR incontinence cover with open back side with height 20 or 30 cm**

90 cm width DKK. 754,-  
140 cm width DKK. 1254,-

**Comfor PUR incontinence sheets**

200 x 140 cm DKK. 526,-  
140 x 100 cm DKK. 279,-

Standard colour is light grey.