

Choose type of cover for SAFE Med Baby mattress and Incubator mattress in the desired size.

Soft-Cell or Comfor PUR Incontinence cover as full cover with zipper or as a cover with open back side.

See product details and prices and state the desired cover type in connection with ordering.

- Soft-Cell™ waterproof cover
 - With fabric surface like SILK
 - Fabric surface with 30g PU waterproof backside
 - Strong, elastic, anti-static, breathable and absorbs more moisture than cotton
 - Ignition test EN 597-1:1995 & ISO 12952-1:2011
 - Washable in machine 60°C-95°C
 - By 60°C – at least 100 times
 - By 95°C – at least 20 times
 - Tumble dry up to 70°C
 - Breathability 1000 g/m² ·24h
 - Öko TEX Class I Product
- Comfor™ PUR Incontinence cover
 - PU waterproof surface on polyester
 - Strong, elastic, breathable
 - Antimicrobial
 - Flametest BS 7175, Crib 5
 - Ignition & match flame test EN 597-1 & 2
 - Washing in machine up to 95°C
 - Tumble dry up to 60°C
 - Test 50 times wash and tumble drying
 - Breathability 500 g/m² ·24h
 - Wipe off with soapy water, alcohol solution or Wet wipe similar to mild chlorine 1000-1200 ppm
 - Biocompatible to ISO 109935

***All covers are FREE from PVC, latex and phthalates**

Prices on Soft-Cell covers are ex. VAT:

Soft-Cell full cover with concealed zipper height 2,5 or 6 cm.

66 x 48 cm	DKK. 339,-
120 x 60 cm	DKK. 805,-
160 x 70 cm	DKK. 845,-
190 x 80 cm	DKK. 868,-

Soft-Cell cover with open back side height 6 cm

70 x 45 cm	DKK. 245,-
120 x 60 cm	DKK. 599,-
160 x 70 cm	DKK. 610,-
190 x 80 cm	DKK. 622,-

Soft-Cell sheet

200 x 140 cm	DKK. 526,-
140 x 100 cm	DKK. 279,-

Standard colour is white.

Prices on Comfor PUR incontinence covers are ex. VAT:

Comfor PUR incontinence full cover with concealed plastic zipper and glued seams height 2,5 cm or 6 cm.

66 x 48 cm	kr. 719,-	sort farve
------------	-----------	------------

Comfor PUR incontinence full cover with concealed plastic zipper height 6 cm, light grey colour

120 x 60 cm	DKK. 805,-
160 x 70 cm	DKK. 845,-
190 x 80 cm	DKK. 868,-

Comfor PUR incontinence cover with open back side in height 6 cm

70 x 45 cm	DKK. 245,-
120 x 60 cm	DKK. 599,-
160 x 70 cm	DKK. 610,-
190 x 80 cm	DKK. 622,-

Comfor sheet

200 x 140 cm	DKK. 526,-
140 x 100 cm	DKK. 279,-

Standard colour on cover for incubator is black.

Standard colour on sheet and cover for baby bed is light grey.