




Comply ignition & match flame test

Cigaret	<b>BS EN 597-1</b>
Match flame	<b>BS EN 597-2</b>
Cigaret	<b>BS EN 1021-1</b>
Match flame	<b>BS EN 1021-2</b>

## SAFE Med OP-positioning cushions and mattresses in SAFE ViscoPore foam and TEX foam with Comfor OP cover



SAFE ViscoPore foam is visco-elastic temperature sensitive foam shaping to the body. The foam is pressure distributing, soundless and provides stability. The foam is firm by low temperatures and by body warmth and pressure it shapes to the body and distributes pressure.

TEX is a stable elastic foam used for specific cushions.

All foams are OEKO-TEX® class 1 product.

Ignition and – match flame test EN 1021-1 & 2, EN 597-1 & 2  
CFC free, halogen free, phthalates free, latex free

Breathable – absorbs and wicks away moisture from the body.

Machine wash test 10 times at 95°C, tumble dryer test 10 times at 90°C, autoclave test 10 times at 134°C

Anti-bacterial, anti-static, 100% allergy friendly


CE-marked medical device MDR-regulation (EU) 2017/745



## COMFOR™ OP Cover - black color

- Cover for cushions & mattress to the OP-theater, 210g/m<sup>2</sup>
- Polyester 125g, PU coating 85 g  
Free from PVC, latex and phthalates  
Biocompatible to ISO 10993-5
- Tested as medical device
- Strong, elastic, breathable
- Antimicrobiel
- Permanently anti-static (BS 2050: ≤ 10<sup>6</sup> Ohm)
- X-ray penetrable and it is harmless by occasional contact
- Flametest BS 7175, Crib 5  
Ignition and match flame test EN 597-1 & 2
- Machine wash up to 95°C / test 50 times
- Tumble dry up to 60°C / test 50 times
- Wipe off with soapy water, alcohol solution or Chlorine  
Disinfection or similar with active chlorine 1000-1200 ppm
- 200 cm water-column(waterproof)
- CE-marked Medical Device MDR regulation (EU) 2017/745

**Comfor OP covers have concealed plastic zipper and glued seams.  
Can be supplied with anti-slip backside**



# SAFE-MED

Pressure Distributing Products **OP**

## Back position

No. 122OP  
Pillow  
for stability  
50x31x9,5/7,5 cm



No. 115OP  
Pillow for total  
pressure relief of  
back head  
40x32x6 cm



No. 66.01OP  
Head fixation  
33x29x12 cm



No. 121OP  
Neck pillow with  
side support  
27x20x10/5/10 cm



No. 77.00OP  
Neck roll pillow  
25x12 cm Ø



No. 11.01OP, 30x40 cm  
No. 11.02OP, 40x50 cm  
No. 11.03OP, 40x60 cm  
Positioning cushion that  
molds to the head and  
neck region



No. 136OP.2  
Arm support Reison  
70x17x10/6,5 cm



No. 136OP.1  
Arm support Nordic  
70x17x10/6,5 cm



No. 128OP, 90x25x6 cm  
No. 129OP, 60x20x6 cm  
No. 129OP.IC, 50x30x2,5 cm  
Leg/Arm support  
- One side with profiled surface  
for pressure distribution



No. 131OP or 131OP.25  
Knee support, 50x19x10 cm  
Knee support, 25x19x10 cm



No. 124OP, 124OP.1 or 124OP.25cm  
Knee support, 50x20x12 cm  
Knee support, 50x20x 8 cm  
Knee support, 25x20x12 cm



No. 133OP  
Cylinder, 50x18 cm Ø



No. 123OP.2  
Leg support-double  
54x20x10/5/10 cm



No. 109OP.1&2  
Leg cushion with  
groove shape for  
heel lift  
30x20x10/7 cm  
60x20x10/7 cm



No. 44.01OP, 160 cm  
No. 44.02OP, 180 cm  
- Long cushion, bended

No. 11.02OP, 40x50 cm  
Positioning cushion





# SAFE-MED

Pressure Distributing Products **OP**

## Side position

No. 122OP  
Pillow for stability  
50x31x9,5/7,5 cm



No. 115OP  
Pillow for total  
pressure relief  
of the ear  
40x32x6 cm



No. 66.01OP  
Head fixation  
33x29x12 cm



No. 106OP.1, 45x30x3/11 cm  
No. 106OP.2, 45x40x3/11 cm  
No. 106OP.3, 80x30x3/12 cm  
No. 106OP.3518 35x18x3/8 cm  
Side wedges for stabilisation



No. 123OP.2 Leg support-double, 54x20x10/5/10 cm  
No. 123OP.1 Leg support-single, 27x20x10/5/10 cm



For positioning of legs in side position.  
18 cm wide x 9 cm high tunnel opening

No. 127OP.25 Tunnel cushion  
30x25x18/9 cm

No. 127OP.50 Tunnel cushion  
30x50x18/9 cm

25 cm wide x 14 cm high tunnel opening

No. 127OP.64 Tunnel cushion  
49x64x21/14 cm

No. 11.02OP,  
Pillow for  
shaping stability  
40x50 cm




No. 11.01OP, 30x40 cm


No. 11.04OP, 40x75 cm

Positioning cushion for shaping stability

No. 125OP  
Side positioning/  
Shoulder support  
85x50x13 cm  
No. 126OP  
Side or supine  
positioning  
cushion including  
insert pad



No. 104OP.2515  
Soft-Lift cushions, set 2 pcs.  
Small cushions to place  
below the patient



## No. 101OP, Mattress for OP bench


\*Can be ordered in other dimensions as required



User weight up to 130 kg

**No.101OP.190, 190x50x6 cm - Yellow/Pink SAFE ViscoPore**  
**No. 101OP.95, 95x50x6 cm**

# Prone Positioning System



2 cm height of base

Prone positioning system 3 parts, no. 135OP  
Width 54 cm and height 20 cm

Support below feet with no. 131OP

Head support no. 134OP  
29 x 24 x 16 cm  
Length x width x height

\*The hole for tubes is right to the bottom.

Knee cushion for pressure relief below the knees  
no. 105YGOP.2020  
20 x 20 x 7 cm





# SAFE-MED

Pressure Distributing Products **OP**

- No. 101ICOP.1015, 10x15x2,5 cm Yellow soft, flat
- No. 101ICOP.3020, 30x20x2,5 cm Yellow soft, flat
- No. TEXOP.3020, 30x20x6 cm TEX, flat
- No. TEX101OP.3020, 30x20x9 cm TEX + Pink medium, flat
- No. 105YGOP.2020 20x20x7 cm Yellow/Grey 3+4 cm soft/firm, profiled
- No. 201YGOP.3050 30x50x10 cm Yellow/Grey 5+5 cm soft/firm, flat



Rectangular positioning cushions in various sizes and constructions that distribute pressure and are stable for legs, arms, knees with Comfor OP black cover. Cover is with 100% tight glued seams and concealed plastic zipper. Can be supplied with anti-slip back side.

\*The cushions can be ordered in other dimensions.

Pictures show some cushions with covers and how some cushions look without covers.

The cushions shape exactly to the body supplying stability, support, comfort and pain relief

Silent material – no noise by movements

No deformation of tissue:

Anti-decubitus - perfect pressure distribution for prevention and healing of pressure sores



## The Pressure Distributing Positioning System

### Type 10, 11 cushions

Rectangular cushions

No. 10.00OP	15 x 15 x 5 cm
No. 11.00OP	18 x 30 cm
No. 11.01OP	30 x 40 cm
No. 11.02OP	40 x 50 cm
No. 11.03OP	40 x 60 cm
No. 11.04OP	40 x 75 cm
Height app. 7/8 cm	

### Type 22 cushions

Rectangular cushions

No. 22.01OP	50 x 60 cm
No. 22.02OP	50 x 90 cm
No. 22.03OP	60 x 90 cm
No. 22.04OP	65 x 100 cm
Height app. 11 cm	

### Type 33 cushion

Boomerang

No. 33.01OP	160 cm
	26 (40) cm
	Height 10 cm
	70°

### Type 44 cushions

Long Cushions

No. 44.01OP	160 x 30 cm
No. 44.02OP	180 x 30 cm
No. 44.04OP	220 x 30 cm

Neck support  
No. 66.01OP  
33 x 29 x 12cm

Neck roll pillow  
No. 77.00OP  
25 cm Ø 12 cm



Filling is 100% viscoelastic foam granulate: SAFE ViscoPore / Fabric material is 100% Cotton  
Cushions can be **machine washed up to 95°C and tumble dried. Optimal hygiene and cleanliness can be maintained.** Cover to wipe off with soapy water, alcohol solution or Chlorine Disinfection or similar with active chlorine 1000-1200 ppm.