

Choose type of cover for Harmony cushions no. 10-11.00,11.01,11.02,11.03,11.04 : Soft-Cell, Comfor PUR or Cotton velour. See product details and prices and state the desired cover type in connection with ordering.

- Soft-Cell™ waterproofcover
- With fabric surface like SILK
- Fabric surface with 30g PU waterproofbackside
- Strong, elastic, anti-static, breathable and absorbs more moisture than cotton
- Ignition test EN 597-1:1995 & ISO 12952-1:2011
- Washable in machine 60°C-95°C
- By 60°C – at least 100 times
- By 95°C – at least 20 times
- Tumble dry up to 70°C
- Breathability 1000 g/m² ·24h
- Öko TEX Class I Product
- Comfor™ PUR Incontinence cover
- PU waterproof surface on polyester
- Strong, elastic, breathable
- Antimicrobiel
- Flametest BS 7175, Crib 5
- Ignition & match flame test EN 597-1 & 2
- Washing in machine up to 95°C
- Tumble dry up to 60°C
- Test 50 times wash and tumble drying
- Breathability 500 g/m² ·24h
- Wipe off with soapy water, alcohol solution or Wet wipe similar to mild chlorine 1000-1200 ppm
- Biocompatible to ISO 109935

***All covers are FREE from PVC, latex and phthalates**

Prices on Soft-Cell covers are ex. VAT:

Soft-Cell cover with concealed zipper

15 x 15 x 5 cm	DKK. 163,-
18 x 30 x 8 cm	DKK. 178,-
30 x 40 x 8 cm	DKK. 191,-
40 x 50 x 8 cm	DKK. 200,-
40 x 60 x 8 cm	DKK. 223,-
40 x 75 x 8 cm	DKK. 266,-

Standard colour on cover is white.

Colour green, blue, grey, black to choose + DKK. 50,-

Prices on Comfor PUR incontinence cover are ex. VAT:

Comfor PUR incontinence cover with concealed zipper

15 x 15 x 5 cm	DKK. 141,-
18 x 30 x 8 cm	DKK. 149,-
30 x 40 x 8 cm	DKK. 174,-
40 x 50 x 8 cm	DKK. 226,-
40 x 60 x 8 cm	DKK. 243,-
40 x 75 x 8 cm	DKK. 293,-

Standard colour on cover is light grey.

Cotton velour cover, very soft and elastic

- 80% Cotton-20% Polyester
- Machine wash 60°C – no shrinkage
- Machine wash 95°C – up to 5% shrinkage
- Tumble dry at medium temperature

Prices on Cotton velour cover are ex. VAT:

Cotton velour cover with overlap closing

15 x 15 x 5 cm	DKK. 141,-
18 x 30 x 8 cm	DKK. 145,-
30 x 40 x 8 cm	DKK. 157,-
40 x 50 x 8 cm	DKK. 197,-
40 x 60 x 8 cm	DKK. 242,-
40 x 75 x 8 cm	DKK. 270,-

Standard colour on cover is white.