

Choose type of cover for SAFE Med seat cushion no. 104, 104YP, 105P, 105G, 105YG in the desired size.

See product details and prices and state the desired cover type in connection with ordering.

- Soft-Cell™ waterproof cover
- With fabric surface like SILK
- Fabric surface with 30g PU waterproof backside
- Strong, elastic, anti-static, breathable and absorbs more moisture than cotton
- Ignition test EN 597-1:1995 & ISO 12952-1:2011
- Washable in machine 60°C-95°C
- By 60°C – at least 100 times
- By 95°C – at least 20 times
- Tumble dry up to 70°C
- Breathability 1000 g/m² · 24h
- Öko TEX Class I Product
- Comfor™ PUR Incontinence cover
- PU waterproof surface on polyester
- Strong, elastic, breathable
- Antimicrobial
- Flametest BS 7175, Crib 5
- Ignition & match flame test EN 597-1 & 2
- Washing in machine up to 95°C
- Tumble dry up to 60°C
- Test 50 times wash and tumble drying
- Breathability 500 g/m² · 24h
- Wipe off with soapy water, alcohol solution or Wet wipe similar to mild chlorine 1000-1200 ppm
- Biocompatible to ISO 109935

***All covers are FREE from PVC, latex and phthalates**

Prices on Soft-Cell covers are ex. VAT:

Soft-Cell full cover with concealed zipper on 2 sides, anti-slip back side and handle

40 x 40 cm	DKK. 305,-
45 x 40 cm	DKK. 305,-
45 x 45 cm	DKK. 317,-
50 x 45 cm	DKK. 317,-
50 x 50 cm	DKK. 317,-
50 x 55 cm	DKK. 328,-
55 x 55 cm	DKK. 328,-
60 x 50 cm	DKK. 339,-
65 x 65 cm	DKK. 339,-

Soft-Cell cover with open back side

40 x 40 cm	DKK. 163,-
45 x 40 cm	DKK. 163,-
45 x 45 cm	DKK. 174,-
50 x 45 cm	DKK. 174,-
50 x 50 cm	DKK. 174,-

Standard colour on covers is black.

Anti-slip back side is grey.

Prices on Comfor PUR incontinence cover are ex. VAT:

Comfor PUR incontinence full cover with concealed zipper on 2 sides, anti-slip back side and handle

40 x 40 cm	DKK. 347,-
45 x 40 cm	DKK. 347,-
45 x 45 cm	DKK. 389,-
50 x 45 cm	DKK. 317,-
50 x 50 cm	DKK. 317,-
50 x 55 cm	DKK. 328,-
55 x 55 cm	DKK. 328,-
60 x 50 cm	DKK. 339,-
65 x 65 cm	DKK. 339,-

Comfor PUR incontinence cover with open back side

40 x 40 cm	DKK. 181,-
45 x 40 cm	DKK. 181,-
45 x 45 cm	DKK. 194,-
50 x 45 cm	DKK. 194,-
50 x 50 cm	DKK. 194,-

Standard colour on covers is black.

Anti-slip back side is grey.