

**Choose type of cover for SAFE Med Wedge no. 106 in the desired size.**

See product details and prices and state the desired cover type in connection with ordering.

- Soft-Cell™ waterproof cover
- With fabric surface like SILK
- Fabric surface with 30g PU waterproof backside
- Strong, elastic, anti-static, breathable and absorbs more moisture than cotton
- Ignition test EN 597-1:1995 & ISO 12952-1:2011
- Washable in machine 60°C-95°C
- By 60°C – at least 100 times
- By 95°C – at least 20 times
- Tumble dry up to 70°C
- Breathability 1000 g/m<sup>2</sup> · 24h
- Öko TEX Class I Product
- Comfor™ PUR Incontinence cover
- PU waterproof surface on polyester
- Strong, elastic, breathable
- Antimicrobial
- Flametest BS 7175, Crib 5
- Ignition & match flame test EN 597-1 & 2
- Washing in machine up to 95°C
- Tumble dry up to 60°C
- Test 50 times wash and tumble drying
- Breathability 500 g/m<sup>2</sup> · 24h
- Wipe off with soapy water, alcohol solution or Wet wipe similar to mild chlorine 1000-1200 ppm
- Biocompatible to ISO 109935

**\*All covers are FREE from PVC, latex and phthalates**

**Prices on Soft-Cell covers are ex. VAT:**

**Soft-Cell cover with concealed zipper on 2 sides and anti-slip back side**

45 x 30 x 3/11 cm	DKK. 333,-
45 x 40 x 3/11 cm	DKK. 333,-
80 x 30 x 3/12 cm	DKK. 596,-

**Soft-Cell cover with open back side**

45 x 30 x 3/11 cm	DKK. 218,-
45 x 40 x 3/11 cm	DKK. 218,-

**Standard colour on cover is white.**

**Anti-slip back side is grey.**

**Prices on Comfor PUR incontinence cover are ex. VAT:**

**Comfor PUR incontinence cover with concealed zipper on 2 sides and anti-slip back side**

45 x 30 x 3/11 cm	DKK. 333,-
45 x 40 x 3/11 cm	DKK. 333,-
80 x 30 x 3/12 cm	DKK. 596,-

**Comfor PUR incontinence cover with open back side**

45 x 30 x 3/11 cm	DKK. 218,-
45 x 40 x 3/11 cm	DKK. 218,-

**Standard colour on cover is light grey.**

**Anti-slip back side is grey.**