

**Choose type of cover for SAFE Med Leg and Arm support no. 109 in the desired size.**

See product details and prices and state the desired cover type in connection with ordering.

- Soft-Cell™ waterproof cover
- With fabric surface like SILK
- Fabric surface with 30g PU waterproof backside
- Strong, elastic, anti-static, breathable and absorbs more moisture than cotton
- Ignition test EN 597-1:1995 & ISO 12952-1:2011
- Washable in machine 60°C-95°C
- By 60°C – at least 100 times
- By 95°C – at least 20 times
- Tumble dry up to 70°C
- Breathability 1000 g/m<sup>2</sup> ·24h
- Öko TEX Class I Product
- Comfor™ PUR Incontinence cover
- PU waterproof surface on polyester
- Strong, elastic, breathable
- Antimicrobial
- Flametest BS 7175, Crib 5
- Ignition & match flame test EN 597-1 & 2
- Washing in machine up to 95°C
- Tumble dry up to 60°C
- Test 50 times wash and tumble drying
- Breathability 500 g/m<sup>2</sup> ·24h
- Wipe off with soapy water, alcohol solution or Wet wipe similar to mild chlorine 1000-1200 ppm
- Biocompatible to ISO 10993-5

**\*All covers are FREE from PVC, latex and phthalates**

**Prices on Soft-Cell covers are ex. VAT:**

**Soft-Cell cover with concealed zipper on 2 sides and anti-slip back side**

30 x 20 x 10/7/10 cm	DKK. 318,-
40 x 20 x 10/7/10 cm	DKK. 325,-
50 x 20 x 10/7/10 cm	DKK. 325,-
60 x 20 x 10/7/10 cm	DKK. 325,-

**Standard colour on cover is white.**

**Anti-slip back side is grey.**

**Prices on Comfor PUR incontinence cover are ex. VAT:**

**Comfor PUR incontinence cover with concealed zipper and anti-slip back side**

30 x 20 x 10/7/10 cm	DKK. 318,-
40 x 20 x 10/7/10 cm	DKK. 325,-
50 x 20 x 10/7/10 cm	DKK. 325,-
60 x 20 x 10/7/10 cm	DKK. 325,-

**Standard colour on cover is light grey.**

**Anti-slip back side is grey.**