

Choose type of cover for SAFE Med Ear pillow no. 115: Soft-Cell, Comfor PUR or Cotton velour

Soft-Cell cover with hole and concealed zipper is the most used cover for ear pillow.

See product details and prices and state the desired cover type in connection with ordering.

- Soft-Cell™ waterproof cover
- With fabric surface like SILK
- Fabric surface with 30g PU waterproof backside
- Strong, elastic, anti-static, breathable and absorbs more moisture than cotton
- Ignition test EN 597-1:1995 & ISO 12952-1:2011
- Washable in machine 60°C-95°C
- By 60°C – at least 100 times
- By 95°C – at least 20 times
- Tumble dry up to 70°C
- Breathability 1000 g/m² ·24h
- Öko TEX Class I Product
- Comfor™ PUR Incontinence cover
- PU waterproof surface on polyester
- Strong, elastic, breathable
- Antimicrobial
- Flametest BS 7175, Crib 5
- Ignition & match flame test EN 597-1 & 2
- Washing in machine up to 95°C
- Tumble dry up to 60°C
- Test 50 times wash and tumble drying
- Breathability 500 g/m² ·24h
- Wipe off with soapy water, alcohol solution or Wet wipe similar to mild chlorine 1000-1200 ppm
- Biocompatible to ISO 10993-5

***All covers are FREE from PVC, latex and phthalates**

Prices on Soft-Cell covers are ex. VAT:

Soft-Cell cover with hole and concealed zipper

40 x 32 x 6 cm DKK. 233,-

Soft-Cell cover without hole and overlap closing

40 x 32 x 6 cm DKK. 107,-

Standard colour on cover is white.

Prices on Comfor PUR incontinence cover are ex. VAT:

Comfor PUR incontinence cover with concealed zipper

40 x 32 x 6 cm DKK. 318,-

Standard colour on cover is light grey.

Cotton velour cover, very soft and elastic

- 80% Cotton-20% Polyester
- Machine wash 60°C – no shrinkage
- Machine wash 95°C – up to 5% shrinkage
- Tumble dry at medium temperature

Prices on Cotton velour cover are ex. VAT:

Cotton velour cover with hole and concealed zipper

40 x 32 x 6 cm DKK. 318,-

Standard colour on cover is white.