

**Choose type of cover for SAFE Med Knee support no. 124, 124.1 & 131 in the desired size.**

See product details and prices and state the desired cover type in connection with ordering.

- Soft-Cell™ waterproofcover
- With fabric surface like SILK
- Fabric surface with 30g PU waterproofbackside
- Strong, elastic, anti-static, breathableand absorbs more moisture than cotton
- Ignition test EN 597-1:1995 & ISO 12952-1:2011
- Washable in machine 60°C-95°C
- By 60°C – at least 100 times
- By 95°C – at least 20 times
- Tumble dry up to 70°C
- Breathability1000 g/m2 ·24h
- Öko TEX Class I Product
- Comfor™ PUR Incontinence cover
- PU waterproofsurface on polyester
- Strong, elastic, breathable
- Antimicrobiel
- Flametest BS 7175, Crib 5
- Ignition & match flame test EN 597-1 & 2
- Washing in machine up to 95°C
- Tumble dry up to 60°C
- Test 50 times wash and tumbledrying
- Breathability500 g/m2 ·24h
- Wipe off with soapy water, alcohol solution orWet wipe similar to mild chlorine 1000-1200 ppm
- Biocompatible to ISO 109935

**\*All covers are FREE from PVC, latex and phthalates**

**Prices on Soft-Cell covers are ex. VAT:**

**Soft-Cell cover with concealed zipper on 2 sides and anti-slip back side**

50 x 20 x 12 cm	DKK. 387,-
50 x 20 x 8 cm	DKK. 387,-
50 x 19 x 10 cm	DKK. 387,-

**Standard colour on cover is white.**

**Black colour can be choosen.**

**Anti-slip back side is grey.**

**Prices on Comfor PUR incontinence cover are ex. VAT:**

**Comfor PUR incontinence cover with concealed zipper and anti-slip back side**

50 x 20 x 12 cm	DKK. 387,-
50 x 20 x 8 cm	DKK. 387,-
50 x 19 x 10 cm	DKK. 387,-

**Standard colour on cover is light grey.**

**Anti-slip back side is grey.**