

Choose type of cover for SAFE Med Leg and Arm rest no. 136.50 & 136.70 in the desired size.

See product details and prices and state the desired cover type in connection with ordering.

- Soft-Cell™ waterproofcover
- With fabric surface like SILK
- Fabric surface with 30g PU waterproofbackside
- Strong, elastic, anti-static, breathableand absorbs more moisture than cotton
- Ignition test EN 597-1:1995 & ISO 12952-1:2011
- Washable in machine 60°C-95°C
- By 60°C – at least 100 times
- By 95°C – at least 20 times
- Tumble dry up to 70°C
- Breathability1000 g/m2 ·24h
- Öko TEX Class I Product
- Comfor™ PUR Incontinence cover
- PU waterproofsurface on polyester
- Strong, elastic, breathable
- Antimicrobiel
- Flametest BS 7175, Crib 5
- Ignition & match flame test EN 597-1 & 2
- Washing in machine up to 95°C
- Tumble dry up to 60°C
- Test 50 times wash and tumbledrying
- Breathability500 g/m2 ·24h
- Wipe off with soapy water, alcohol solution orWet wipe similar tomild chlorine 1000-1200 ppm
- Biocompatible to ISO 109935

***All covers are FREE from PVC, latex and phthalates**

Prices on Soft-Cell covers are ex. VAT:

Soft-Cell cover with concealed zipper on 2 sides and anti-slip back side

50 x 17,5 x 10/6,5/10 cm DKK. 325,-
70 x 17,5 x 10/6,5/10 cm DKK. 325,-

Standard colour on cover is white.

Anti-slip back side is grey.

Prices on Comfor PUR incontinence cover are ex. VAT:

Comfor PUR incontinence cover with concealed zipper and anti-slip back side

50 x 17,5 x 10/6,5/10 cm DKK. 325,-
70 x 17,5 x 10/6,5/10 cm DKK. 325,-

Standard colour on cover is light grey.

Anti-slip back side is grey.