

Choose type of cover for Harmony cushions no. 22.01, 22.02,22.03, 22.04: Soft-Cell, Comfor PUR or Cotton velour.

See product details and prices and state the desired cover type in connection with ordering.

- Soft-Cell™ waterproofcover
- With fabric surface like SILK
- Fabric surface with 30g PU waterproofbackside
- Strong, elastic, anti-static, breathable and absorbs more moisture than cotton
- Ignition test EN 597-1:1995 & ISO 12952-1:2011
- Washable in machine 60°C-95°C
- By 60°C – at least 100 times
- By 95°C – at least 20 times
- Tumble dry up to 70°C
- Breathability 1000 g/m² ·24h
- Öko TEX Class I Product
- Comfor™ PUR Incontinence cover
- PU waterproof surface on polyester
- Strong, elastic, breathable
- Antimicrobiel
- Flametest BS 7175, Crib 5
- Ignition & match flame test EN 597-1 & 2
- Washing in machine up to 95°C
- Tumble dry up to 60°C
- Test 50 times wash and tumble drying
- Breathability 500 g/m² ·24h
- Wipe off with soapy water, alcohol solution or Wet wipe similar to mild chlorine 1000-1200 ppm
- Biocompatible to ISO 109935

***All covers are FREE from PVC, latex and phthalates**

Prices on Soft-Cell covers are ex. VAT:

Soft-Cell cover with concealed zipper

50 x 60 x 11 cm	DKK. 264,-
50 x 90 x 11 cm	DKK. 307,-
60 x 90 x 11 cm	DKK. 358,-
65 x 100 x 11 cm	DKK. 421,-

Standard colour on cover is white.

Colour green, blue, grey, black to choose + DKK. 50,-

Prices on Comfor PUR incontinence cover are ex. VAT:

Comfor PUR incontinence cover with concealed zipper

50 x 60 x 11 cm	DKK. 290,-
50 x 90 x 11 cm	DKK. 368,-
60 x 90 x 11 cm	DKK. 458,-
65 x 100 x 11 cm	DKK. 475,-

Standard colour on cover is light grey.

Cotton velour cover, very soft and elastic

- 80% Cotton-20% Polyester
- Machine wash 60°C – no shrinkage
- Machine wash 95°C – up to 5% shrinkage

Tumble dry at medium temperature

Prices on Cotton velour cover are ex. VAT:

Cotton velour cover with overlap closing

50 x 60 x 11 cm	DKK. 268,-
50 x 90 x 11 cm	DKK. 321,-
60 x 90 x 11 cm	DKK. 400,-
65 x 100 x 11 cm	DKK. 417,-

Standard colour on cover is white.