

**Choose type of cover for Harmony Long cushions no. 44.01, 44.02, 44.04 : Soft-Cell, Comfor PUR or Cotton velour.**

See product details and prices and state the desired cover type in connection with ordering.

- Soft-Cell™ waterproofcover
- With fabric surface like SILK
- Fabric surface with 30g PU waterproofbackside
- Strong, elastic, anti-static, breathable and absorbs more moisture than cotton
- Ignition test EN 597-1:1995 & ISO 12952-1:2011
- Washable in machine 60°C-95°C
- By 60°C – at least 100 times
- By 95°C – at least 20 times
- Tumble dry up to 70°C
- Breathability 1000 g/m<sup>2</sup> ·24h
- Öko TEX Class I Product
- Comfor™ PUR Incontinence cover
- PU waterproof surface on polyester
- Strong, elastic, breathable
- Antimicrobiel
- Flametest BS 7175, Crib 5
- Ignition & match flame test EN 597-1 & 2
- Washing in machine up to 95°C
- Tumble dry up to 60°C
- Test 50 times wash and tumble drying
- Breathability 500 g/m<sup>2</sup> ·24h
- Wipe off with soapy water, alcohol solution or Wet wipe similar to mild chlorine 1000-1200 ppm
- Biocompatible to ISO 109935

**\*All covers are FREE from PVC, latex and phthalates**

**Prices on Soft-Cell covers are ex. VAT:**

**Soft-Cell cover with concealed zipper**

160 x 30 x 15 cm	DKK. 457,-
180 x 30 x 15 cm	DKK. 478,-
220 x 30 x 15 cm	DKK. 526,-

**Standard colour on cover is white.**

**Colour green, blue, grey, black to choose + DKK. 50,-**

**Prices on Comfor PUR incontinence cover are ex. VAT:**

**Comfor PUR incontinence cover with concealed zipper**

160 x 30 x 15 cm	DKK. 498,-
180 x 30 x 15 cm	DKK. 500,-
220 x 30 x 15 cm	DKK. 544,-

**Standard colour on cover is light grey.**

**Cotton velour cover, very soft and elastic**

- 80% Cotton-20% Polyester
- Machine wash 60°C – no shrinkage
- Machine wash 95°C – up to 5% shrinkage

**Tumble dry at medium temperature**

**Prices on Cotton velour cover are ex. VAT:**

**Cotton velour cover with overlap closing**

160 x 30 x 15 cm	DKK. 437,-
180 x 30 x 15 cm	DKK. 481,-
220 x 30 x 15 cm	DKK. 526,-

**Standard colour on cover is white.**